

SNS-BC




SAFETY WINGED BLOOD COLLECTION SET

**WORLD-CLASS SAFETY WINGED BLOOD COLLECTION DEVICE
WITH OR WITHOUT LUER ADAPTER**



This world-class safety winged blood collection device offers a painless venipuncture and optimal blood flow, ensuring safety for the patient and healthcare professionals.



ULTRA THIN WALLED NEEDLE



FLEXIBLE WINGS BUT RIGID ENOUGH AND TEXTURED TO ENSURE A GOOD GRIP

INTUITIVE USE



IN VEIN ACTIVATION WITH ONE HAND



TUBING WITH NO MEMORY



BILINGUAL AND EASY TO USE PACKAGING



STERILE/LATEX FREE/DEHP FREE

INSTRUCTIONS FOR USE & REORDER NUMBERS



Step1:
Venipuncture.



Step2:
Insert the tube for
blood collection.



Step3:
After completion of blood
collection, use middle finger to
cover the puncture site with a
gauze pad and place the index
finger on the curved design to
fix, the other hand to withdraw
the tubing back to activate the
safety mechanism.



Step4:
Audible click indicates
activation.
Disposal after use.

DESCRIPTION	REORDER NUMBER	IMAGE
SNS BC 21G 3/4" needle 10cm tubing	SNS-LC1021	
SNS BC 23G 3/4" needle 10cm tubing	SNS-LC1023	
SNS BC 25G 3/4" needle 10cm tubing	SNS-LC1025	
SNS BC 21G 3/4" needle 19cm tubing	SNS-LC1921	
SNS BC 23G 3/4" needle 19cm tubing	SNS-LC1923	
SNS BC 25G 3/4" needle 19cm tubing	SNS-LC1925	
SNS BC 21G 3/4" needle 30cm tubing	SNS-LC3021	
SNS BC 23G 3/4" needle 30cm tubing	SNS-LC3023	
SNS BC 25G 3/4" needle 30cm tubing	SNS-LC3025	
SNS BC with Luer adapter 21G 3/4" needle 19cm tubing	SNS-LA1921	
SNS BC with Luer adapter 23G 3/4" needle 19cm tubing	SNS-LA1923	
SNS BC with Luer adapter 25G 3/4" needle 19cm tubing	SNS-LA1925	
SNS BC with Luer adapter 21G 3/4" needle 30cm tubing	SNS-LA3021	
SNS BC with Luer adapter 23G 3/4" needle 30cm tubing	SNS-LA3023	
SNS BC with Luer adapter 25G 3/4" needle 30cm tubing	SNS-LA3025	
SNS BC with holder 21G 3/4" needle 19cm tubing	SNS-LH1921	
SNS BC with holder 23G 3/4" needle 19cm tubing	SNS-LH1923	
SNS BC with holder 25G 3/4" needle 19cm tubing	SNS-LH1925	
SNS BC with holder 21G 3/4" needle 30cm tubing	SNS-LH3021	
SNS BC with holder 23G 3/4" needle 30cm tubing	SNS-LH3023	
SNS BC with holder 25G 3/4" needle 30cm tubing	SNS-LH3025	